新冠抗原检测试剂产品手册
**Product Manual of
COVID-19 Antigen Detection Kit**

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COMPANY PROFILE

New Gene (Hangzhou) Bioengineering Co., Ltd. is located in Hangzhou, China. It is a high-tech company engaged in the research, development, manufacture and distribution of biological products. It is committed to creating biological materials such as antigens and antibodies, in vitro diagnostic reagents and related devices, and also the complete industrial chain of artificial intelligence assisted diagnosis system. The product line covers a full range of in vitro diagnostic products such as immune diagnosis, molecular diagnosis, and microbiological testing. NEWGENE has profound technical accumulation and unique technological advantages in the areas of early cancer screening, rapid detection of infectious diseases, and rapid screening of geriatric diseases.




NEWGENE's manufacturing system meets GMP standards for medical devices, and is certified with ISO13485 by British BSI. Relevant in vitro diagnostic reagent products have obtained the EU CE certification. NEWGENE is also a member on the "allow list" issued by Chinese Ministry of Commerce for anti-epidemic products exporting. At present, NEWGENE COVID-19 Antigen Detection Kit has **registered in** many countries, including **Germany, France, Italy, Switzerland, Belgium, Portugal, Czech Republic, Hungary, Greece, Moldova, Peru, Argentina, Ecuador, Kenya, Zimbabwe etc.** NEWGENE COVID-19 Antigen Detection Kit has passed the clinical **validation in** national lab in **Germany, Switzerland, Ecuador, Zimbabwe etc.** NEWGENE products show good performance in sensitivity and specificity compared with international brand products. The products have exported to more than 50 countries and regions.

COVID-19 Antigen Detection Kit

Nasopharyngeal Swab Sample

25 Tests/Box

1 Test/Box


Product Feature

- Samples: **Nasopharyngeal Swab**
- Multiple Packaging Specifications: 25 Tests/Box or 1 Test/Box
- Fast Detection: Result in 15 minutes.
- High Accuracy.
- Easy to Use.


Components

NO.	Components name	25 Tests/Box	1Test/Box
1	Test Card	*25	*1
2	Sample Extraction Tube	*25	*1
3	Tube Cap	*25	*1
4	Sampling Swab	*25	*1
5	Package Insert	*1	*1


SAMPLING STEPS

DETECTION STEPS

INTERPRETATION OF RESULTS

Positive(+): Red bands appear at both of T and C line in 15 to 30 minutes. A white band at the T line should be considered as a negative result.


Negative(-): A red band appears at C line while no red band appears at T line in 15 to 30 minutes after sample loading.

Invalid: As long as no red band appears at C line, it indicates that the test result is invalid, and should retest the sample with another test card.

Sensitivity	Specificity
98.0%	99.1%

Nasal Swab Sample

25 Tests/Box

5 Tests/Box

1 Test/Box


Product Feature

- Samples: **Nasal Swab.**
- Multiple Packaging Specifications: 25 Tests/Box or 1 Test/Box
- Fast Detection: Result in 15 minutes.
- High Accuracy.
- Easy to Use.


Components

NO.	Components name	25 Tests/Box	1Test/Box
1	Test Card	*25	*1
2	Sample Extraction Tube	*25	*1
3	Tube Cap	*25	*1
4	Sampling Swab	*25	*1
5	Package Insert	*1	*1


SAMPLING STEPS

DETECTION STEPS

INTERPRETATION OF RESULTS

Positive(+): Red bands appear at both of T and C line in 15 to 30 minutes. A white band at the T line should be considered as a negative result.

Negative(-): A red band appears at C line while no red band appears at T line in 15 to 30 minutes after sample loading.

Invalid: As long as no red band appears at C line, it indicates that the test result is invalid, and should retest the sample with another test card.

Sensitivity	Specificity
97.1%	99.2%

Oropharyngeal Swab Sample

25 Tests/Box

1 Test/Box


Product Feature

- Samples: **Oropharyngeal Swab** .
- Multiple Packaging Specifications: **25 Tests/Box** or **1 Test/Box**
- Fast Detection: Result in 15 minutes.
- High Accuracy.
- Easy to Use.


Components

NO.	Components name	25 Tests/Box	1Test/Box
1	Test Card	*25	*1
2	Sample Extraction Tube	*25	*1
3	Tube Cap	*25	*1
4	Sampling Swab	*25	*1
5	Package Insert	*1	*1


SAMPLING STEPS

DETECTION STEPS

INTERPRETATION OF RESULTS

Positive(+): Red bands appear at both of T and C line in 15 to 30 minutes. A white band at the T line should be considered as a negative result.

Negative(-): A red band appears at C line while no red band appears at T line in 15 to 30 minutes after sample loading.

Invalid: As long as no red band appears at C line, it indicates that the test result is invalid, and should retest the sample with another test card.

Sensitivity	Specificity
95.7%	99.0%

Sputum Sample (Saliva)

25 Tests/Box

1 Test/Box


Product Feature

- Samples: **Sputum (Saliva).**
- Multiple Packaging Specifications: 25 Tests/Box or 1 Test/Box
- Fast Detection: Result in 15 minutes.
- High Accuracy.
- Easy to Use.


Components

NO.	Components name	25 Tests/Box	1Test/Box
1	Test Card	*25	*1
2	Sample Extraction Tube	*25	*1
3	Tube Cap	*25	*1
4	Paper Cup	*25	*1
5	Sputum Dropper	*25	*1


SAMPLING STEPS

DETECTION STEPS

INTERPRETATION OF RESULTS

Positive(+): Red bands appear at both of T and C line in 15 to 30 minutes. A white band at the T line should be considered as a negative result.

Negative(-): A red band appears at C line while no red band appears at T line in 15 to 30 minutes after sample loading.

Invalid: As long as no red band appears at C line, it indicates that the test result is invalid, and should retest the sample with another test card.

Sensitivity	Specificity
97.3%	99.2%

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
A handwritten signature in blue ink, appearing to read 'M.J. van de Velde'.

Dr. M.J. van de Velde

DECLARATION OF CONFORMITY

Regarding In Vitro Diagnostic Directive (98/79/EC)

Manufacturer: New Gene (Hangzhou) Bioengineering Co., Ltd.

Address: Room 1606, Floor 16, Building 5, 688 Bin'an Road, Changhe Street, Binjiang District, Hangzhou City, Zhejiang Province, P. R. China

EC Representative: SUNGO Europe B.V.

Address: Olympisch Stadion 24, 1076DE Amsterdam, Netherlands

Product Name: COVID-19 Antigen Detection Kit

Product Code: COVID-19-NG08

Specification: 25Tests/Box 1Test/Box

Classification: Others (IVDD)

Conformity Assessment

Procedure: Annex III of In Vitro Diagnostic Directive (98/79/EC)

We herewith declare that the above-mentioned products meet the requirements of In Vitro Diagnostic Directive (98/79/EC) and the following harmonized standards.

EN 23640:2015 EN 13640:2002

EN 980:2016 EN 13641:2002

EN ISO 14971:2019 EN ISO 18113-1:2011

EN 13612:2002 EN ISO 18113-4:2011

On behalf of SUNGO Europe office, I confirmed we are EU REP of the company who issue this document.

Signature:

Name/ Position: Mingfu Li / General Manager

Date: 29/09/2020

Place: Hangzhou, Zhejiang, China

Authorized Signature (S)

By Royal Charter

Certificate of Registration

QUALITY MANAGEMENT SYSTEM - ISO 13485:2016

This is to certify that: New Gene (Hangzhou) Bioengineering Co., Ltd.
Room 1606, 16th Floor, No.5 Building
688 Bin'an Road
Binjiang District
Hangzhou
Zhejiang
310052
China

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Holds Certificate No: **MD 729179**

and operates a Quality Management System which complies with the requirements of ISO 13485:2016 for the following scope:

Design and Development, Manufacture and Distribution of In-vitro Diagnostic Rapid Test Kit of Drug Abuse, Manufacture and Distribution of In-vitro Diagnostic Rapid Test Kit of Infectious Diseases.

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For and on behalf of BSI:

Gary E Slack, Senior Vice President - Medical Devices

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Paul-Ehrlich-Institut 

26.03.2021

Comparative evaluation of the sensitivities of SARS-CoV-2 antigen rapid tests

Aim

Comparison of different antigen rapid tests with using identical sample material

Material

Pools from nasopharyngeal and oropharyngeal swabs.

Dry swabs were included in PBS; moist swabs were already included in the transport media of various compositions. Pools are random mixtures obtained from up to 10 samples of comparable CT values diluted 1:10 in negative samples in PBS. The CT values of a pool were determined by means of different PCR assays, and the putative number of RNA copies calculated with the aid of the INSTAND standards. In the case of the PCRs used, a CT value of 25 corresponds to around 10^6 RNA copies/mL. 18 samples each were analysed with CT<25, 23 samples with CT between 25 and 30, and 9 samples with CT>30. The replication of the virus in cell culture was determined as a possible correlate for infectiousness as another characteristic of the samples.

Method

The pools were aliquoted, frozen, shipped, and thawed for evaluation of the tests. For each test, 50 µL of the pool were analysed using the components of the test provided, e.g. swabs. Laboratories participating in the comparative evaluation included the Robert Koch-Institut, the Paul-Ehrlich-Institut, the reference laboratory for coronaviruses (Charité), and the Institute for Microbiology of the German Army (Bundeswehr).

Summary

This comparative evaluation of a large number of SARS-CoV-2 rapid antigen tests (point of care tests; POCT) of different designs and manufacturers with the same sample set allows an overview of the current state of art regarding sensitivity. The results do not allow any conclusions regarding specificity of the tests.

Those POCTs which have up to now been included in the evaluation and have been assessed as reflecting the current state of the art are listed in the table below. Other tests, which were assessed as not reflecting the state of the art were deleted from the list of the Federal Institute for Drugs and Medical Devices (Bundesinstitut für Arzneimittel und Medizinprodukte, BfArM). This comparative evaluation is constantly continued, and the table is amended accordingly.

You should be aware that this comparative evaluation can only cover a random sample of the SARS-CoV-2 rapid antigen tests listed by the BfArM, thus eligible for refunding, and that few other products could not (yet) be taken into account, despite the interests on the part of the manufacturers/distributors.

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Sars-CoV-2 Antigen Test

with

COVID-19 Antigen Detection Kit (Newgene(Hangzhou) Bioengineering)

Department of microbiology LMZ Dr Risch Pregassona

Validation Interval: 01/21 – 02/21

5 Conclusions

The comparable results with the R^{***} Test makes the newgene test suitable to be introduced in Switzerland. It is important to note the better sensitivity of that assay in comparison to R^{***}.

Zimbabwe Clinical Evaluation Report

NATIONAL MICROBIOLOGY REFERENCE LABORATORY

NEWGENECOVID 19 RAPID ANTIGEN TEST EVALUATION AGAINST COVID 19 RT-PCR COMPARISON REPORT

Conclusion

There was 100 % positive percent agreement between the NEWGENE SARS-CoV- 2 antigen test and SARS CoV-2 PCR test done using stored samples at NMRL. NEWGENE kit is recommended for use in Zimbabwe.

Evaluation team

Stanford Mupandasekwa Lab Scientist Signature P.P. D Date 26/2/2021
Boniface Muzividzi Lab Scientist Signature P.P. D Date 26/2/2021

Reviewed by:

Lucia Sisya Quality officer Signature P.P. D Date 26/02/2021
Agnes Juru Chief lab Scientist Signature P.P. D Date 26/2/21

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Registered Countries or White list (Partially)

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Platnost povolení:.

Belgium

https://www.famhp.be/en/human_use/health_products/medical_devices_accessories/covid_19/tests

The screenshot shows the FAMHP website interface. At the top, there are navigation links for 'About the FAMHP', 'Jobs', 'Publications', 'Press', 'Contact', 'Complaints', and 'Webportal'. A search bar is also present. The main navigation menu includes 'Human use', 'Veterinary use', 'Information for the public', and 'Information for professionals'. The 'Human use' section is active, and the breadcrumb trail reads: Home > Human use > Health Products > Medical devices and their accessories > COVID-19 > Tests.

The 'Tests' section contains a list of links:

- New validation procedure for serological tests and antigenic tests
- How should tests be made available in Belgium?
- List of the recommended tests** (highlighted with a red box)
- Questions and answers about at-home antibody tests to check for antibodies against the novel coronavirus (SARS-CoV-2)
- Swabs: compliance verification

Other features on the page include 'Notification of adverse reactions or incidents', 'PIL and SPC of a medicine' (with a search button), and 'News' (with dates 12/02/2021 and 10/02/2021).

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A red arrow points to the 'Code' column with the text: 'Product code, is the unique 13 digit code that a product acquires after its registration into the EOF's database'.

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MINISTRY OF HEALTH

PHARMACY AND POISONS BOARD

(Section 3B(2)(e) of the Pharmacy and Poisons Act, Cap 244 Laws of Kenya)

IN-VITRO DIAGNOSTIC EMERGENCY USE AUTHORIZATION

This Emergency Use Authorization is issued to **New Gene (Hangzhou) Bioengineering Co., Ltd**, for distribution and sale of **Novel Coronavirus Antigen Detection Kit (Colloidal Gold)**

Emergency use Authorization (EUA) No.	MD/2021/7674
EUA valid until	End of COVID -19 Pandemic or EUA revocation
Device category	Medical Device class C/D
GMDN	N/A
GMDN Term	N/A
Intended purpose	For epidemiological COVID-19 Screening
Conditional Approval	N/A

EUA No.: MD/2021/7674

Date of Authorization: March 4, 2021

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Razón Social: **ALCAT S.A.**