

REHAB MOTORIZED TREADMILL

GT8600 RF

Specific designed for disabled and non-disabled people alike. With the increasing recognition that fitness equipment needs to be more accessible and initiative.

Our goal is to ensure that all users have unfettered access to a wide range of exercise environment. With accessibility of GT8600 RF Rehab Treadmill is no longer a barrier to improved health and independence.



USER ADVANTAGES



EN957



8X16 DOT-MATRIX LED CONSOLE

- Reverse Function
- Dual Emergency Stop Buttons and Safety Lanyard
- Heart Rate Monitor
- Intuitive Display
- Visible Iconography
- New-Added 12 Training Programs

FEATURES



A. Dual-Emergency Stop System

To ensure your safety, the treadmill features a double emergency stop system consisting of a safety lanyard and dual emergency stop buttons.

B. Low Step-on Height

Features a low, 16.5 cm/ 6.25" step-on height for easy accessibility.

C. Elevation Level from 0 to 15%

Either forward or backward exercise mode up to 15% power elevation to work different muscle groups.

D. Wide Entry and Accessible for Wheelchair

Consideration of accessibility requirements the running area entrance is with much wider for the limited mobility users to easy get into the treadmill for running.

E. Bi-Directional Belt

Walking or running backward on the treadmill not only will you tone different muscles, but also work on your balance.

F. Wired Controller for Physiotherapist

Good idea to lower the therapy risk during exercising. This unique design is to help physiology trainer to stop the fitness equipment when the user can't be ability for exercise.